



Application Guide



Dermisplus Prevent Pad

With versatile application, the Pads can be used to provide pressure redistribution to many areas of the body including the occiput, between the knees, under the chin or folded for hand contractures.



Dermisplus Prevent Strips

Ideal for protecting pressure points on digits, with medical devices* such as face masks or longer body areas such as the spine.

*Dermisplus Prevent can also be used under medical devices to prevent pressure damage. Example: oxygen masks, nasal specs, tracheostomy tubing and ties, rigid collars.




Dermisplus Prevent Sacrum

Designed for the sacral area, this dressing can also be used around other difficult to dress areas such as ankles or shoulders.



Dermisplus Prevent Heel

Available in two sizes (standard and large), the heel pad can be used long ways across the foot to protect the sole or across the heel to protect the back of the foot. It can also be used to protect elbows.



Effective



Easy Application



Patient Friendly



Cut to size



Washable



Durable

Ordering Details



Description	Size (L/W/H)	Packaging	Order Code
Pad	4 x 4 x 0.1in 10 x 10 x 0.3cm	5/Box	8500007
Pad	4 x 4 x 0.5in 10 x 10 x 1.2cm	5/Box	8500003
Pad	8 x 8 x 0.1in 20 x 20 x 0.3cm	2/Box	8500005
Pad	8 x 8 x 0.5in 20 x 20 x 1.2cm	2/Box	8500009
Strip	12 x 2 x 0.1in 30 x 5 x 0.3cm	5/Box	8500006
Strip	20 x 1 x 0.1in 50 x 2.5 x 0.3cm	5/Box	8500001
Sacrum	10 x 7 x 0.4in 25 x 18.5 x 1cm	1/Box	8500008
Heel	5.7 x 3.7 x 1.6in 14.5 x 9.5 x 4cm	2/Box	8500002
Large Heel	6.4 x 4.7 x 2.7in 16.5 x 12 x 6.8cm	2/Box	8500004

For more information:

Call +1 (716) 713-0549

Email: contactus@frontier-therapeutics.com

Website: www.frontier-therapeutics.com



300 International Drive, Suite 100 Buffalo,
NY 14221 United States of America

© Copyright Frontier Therapeutics Limited, 2019
All rights reserved. GP0132_A70_02

dermisplus[®]
PREVENT